


Auris Hearing Loops: An Overview



Auris Hearing Loops maximize the benefit of hearing aids already fitted (or retro fitted) with “TeleCoils”. The T-Coil receives the wirelessly transmitted signal from the hearing loop: most hearing aids (70% +) are already fitted with T-Coils.

The sound source (mic, TV, soundtrack, mp3 etc) is wirelessly transmitted directly to the recipient’s hearing aid, thus eliminating the problems of:

- 1 Background noise
 - 2 Less than ideal building acoustics
 - 3 Being too far from the sound source
- Too much or too little volume.